

FIZE kUO™

Making Every Drop Count

Urine output measurement and monitoring are vital aspects of critical care and key to early detection of **shock and Acute Kidney Injury** in ICU patients. So how is it possible that we still manage Urine Output manually.



Introducing FIZE kUO™

The technology that puts UO in line with other vital signs



Accurate, minute-to-minute data in real-time from the moment of catheterization



Updates the patient's medical records (EMR) at every given moment



Alerts on changes in patient UO according to personalized thresholds



Frees valuable time of medical staff.



Promotes a hygienic, hands-free environment



Flexible positioning

On or near the bed, not sensitive to tilt



Active drainage mechanism

With auto drainage of tubing to prevent obstruction of flow



Plug and play

Intuitive and easy to operate




Includes a complete closed drainage kit

And can be connected to any size and type of Foley catheter



Suitable for in hospital transfer

4-hour back up battery
MR safe disposable kit



✓ Patented technology

✓ FDA listed

✓ Clinically tested

About FIZEMEDICAL

Founded in 2019, FIZE Medical is the developer of an advanced fluid management platform that brings UO management into the digital healthcare management space. We are passionate about making a difference, and revolutionizing critical care, one drop at a time.

For more information
www.fizemedical.com
info@fizemedical.com